



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

ICU Mechanical Bed ABS Panel and Railing

Product Code :

SAMCOMEDICAL10032A001

Description :

ICU Mechanical Bed ABS Panel and Railing

Technical Specification :

- *Frame work made of rectangular M.S. tube.
- *4 Section top made of perforated M.S. sheet.
- *Back rest, knee rest & TB/RTB positions maneuvered by separate screw from foot end.
- *ABS moulded head and foot panels.
- *ABS four safety side railing.
- * S.S. IV rod with four locations.
- *Bed mounted on 150mm dia castor (2 with brakes).
- *Finish: Epoxy Powder Coated.

Size/Dimension 2120mm X 1000mm X 460-720mm



Samco Medical

Website: www.samcomedical.com,