



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Emergence Vehicle Ambulance

**Product Code :**  
SAMCOMEDICAL10282A001

**Description :**

Emergence Vehicle Ambulance

**Technical Specification :**

Control plate of siren and alarm lamp mounted in driver's cab  
Adjustable volume alarm apparatus  
Red flashing alarm light bar 100w mounted on the front roof of vehicle and with 3 pieces lights  
Red flashing alarm lamps four sides of vehicle outside roof  
Composite material interior decoration in patient compartment  
Medicine cabinet  
Oxygen cylinders 10L  
Ultraviolet radiation disinfector light  
Dirt bucket  
Individual seat  
Long bench one piece (suitable for 3 people) with safe belts, mounted in right side of patient's compartment  
A foldable doctor seat (with safe belt) mounted to floor near to compartment wall  
Partition wall between driver's cab and patient's compartment, a glass window in wall to make observation  
12V light lamp  
Intelligent inverted power supply system (12-220V) 500W  
220V wall outlet  
Interior medical floor  
Cross mark outside marking  
The original air conditioning is installed in driver's cab and in patient's compartment  
Exhaust system is installed on patient's compartment roof  
Multi-functional interior roof in patient's compartment, with handrail, transfusion hanging bracket  
Portable and foldable stretcher  
Automatic loading-and-unloading stretcher  
Hand-free intercom system between driver's cab and patient's compartment  
100W loud speaker  
Oxygen reducing valve  
Fire extinguisher



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),