



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Computer Cardiopulmonary Resuscitation Simulator

Product Code :
SAMCOMEDICAL250001A001

Description :

Computer Cardiopulmonary Resuscitation Simulator

Technical Specification :

Features:

2010 version of the large-screen LCD color display advanced computer cardiopulmonary resuscitation simulator
Large LCD color display:

Artificial respiration and chest compressions, simulated heart beat display, simulated ECG display, huge graph data statistics

Simulate standard airway opening;

Manual chest compressions:

The dynamic bar code indicator shows the pressing depth: the pressing depth is correct (5-6cm area) by the bar code green light, the pressing depth is not enough (less than 5cm) by the bar code yellow, the pressing depth is too deep (greater than 6cm) by the bar code red indicator moving Dynamic feedback shows the depth of CPR compressions.

LCD count display; record the specific reasons for pressing errors in detail (excessive pressing force, too small pressing force, wrong pressing position and correct number of times).

Huge chart data statistics: correct pressing, excessive pressing force, insufficient pressing force, correct pressing, etc. Detailed chart statistics display 4. Language prompt: prompt, detailed prompting of the specific reasons for the pressing error, so that the trainer can correct it in time.

During artificial mouth-to-mouth breathing (breathing):

The dynamic bar code indicator shows the tidal volume: the correct tidal volume (500ml~600ml-1000ml) is displayed by the bar code green light, and the blown tidal volume is too small or too large, respectively, which is moved by the bar code yellow or bar code red indicator light. Feedback shows tidal measurement;

Samco Medical

Website: www.samcomedical.com,

