



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Advanced Endotracheal Intubation Training Model (with Alarm)

**Product Code :**

SAMCOMEDICAL250001A006

**Description :**

Advanced Endotracheal Intubation Training Model (with Alarm)

**Technical Specification :**

Advanced endotracheal intubation training model (with alarm) adopts American Heart Association (AHA) 2005 Cardiopulmonary Resuscitation (CPR) & Cardiovascular Emergency (ECC) guidelines. Under the advice and guidance of domestic authoritative experts, it integrates the advantages of similar foreign products And the joint efforts of the company's scientific and technological personnel, the first to develop a high-end seventh-generation new product specially designed for the popularization of social cardiopulmonary resuscitation and training in medical schools and medical and health systems. Its shape is more perfect, more powerful, and more convenient to operate. .

Main functions advanced endotracheal intubation training model (with alarm):

The function of combining standard human anatomy with intuitive demonstration of real operation.

When performing oral and nasal tracheal intubation training operations, insert the airway correctly and have a side view function; supply air to inflate both lungs, and inject air into the tube balloon to fix the tube.

When the training operation of oral and nasal endotracheal intubation is performed, the wrong operation is inserted into the esophagus, and there is a side visual function and an alarm function. The gas supply distends the stomach.

When the training operation of oral and nasal tracheal intubation is performed, the wrong operation will cause the laryngoscope to cause the teeth to be compressed, and there is an electronic alarm function.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),

