



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Electronic Human Tracheal Intubation Training Model (with Tooth Pressure Alarm)

Product Code :

SAMCOMEDICAL250001A011

Description :

Electronic Human Tracheal Intubation Training Model (with Tooth Pressure Alarm)

Technical Specification :

It can carry out the training operation and teaching demonstration of oral and nasal endotracheal intubation. In the training operation of oral and nasal endotracheal intubation, insert the airway correctly, with electronic display and music function; air supply to inflate both lungs, and inject air into the tube airbag to fix the tube. When the training operation of oral nasal cavity tracheal intubation is performed, the wrong operation is inserted into the esophagus, electronic display and alarm function. The gas supply distends the stomach. When the training operation of oral and nasal tracheal intubation is performed, the wrong operation will cause the laryngoscope to cause the teeth to be compressed, and there is an electronic alarm function. Observe and compare normal pupils on one side and dilated pupils on the other side. Indicates the cricothyroid membrane puncture site.



Samco Medical

Website: www.samcomedical.com,