



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name : Advanced Cardiopulmonary Resuscitation, AED Defibrillation Simulator	Product Code : SAMCOMEDICAL250001A019
--------------------------------------------------------------------------------------------------	-------------------------------------------------

Description :

Advanced Cardiopulmonary Resuscitation, AED Defibrillation Simulator (computer control, two-in-one combination)

Technical Specification :

Heart Association (AHA) 2010 Cardiopulmonary Resuscitation (CPR) & Emergency Cardiovascular (ECC) Guidelines

Main functions of the system: (The simulator can use the existing resource office computer to install the software during the training time. It will not affect other daily operations.)

Cardiopulmonary resuscitation: Advanced cardiopulmonary resuscitation, AED defibrillation simulator (computer control, two-in-one combination)

Simulation of vital signs: pupillary changes and carotid artery pulsation.

Simulates standard airway opening.

Cardiopulmonary resuscitation training, designed according to the 2010 Cardiopulmonary Resuscitation Guidelines, can perform artificial respiration and chest compressions, and prompts throughout the process.

Standard airway opening;

Language prompts for correct and incorrect operations;

Statistical data print results;

You can choose training and popular assessment, professional assessment methods.

Multimedia teaching can be carried out, and it can also be connected to a projector to play teaching.

Student management: objective evaluation report recording and playback, practice and assessment.

AED defibrillation function:

The simulated external defibrillator (special for training) is used, patch electrodes are provided, and prompts are provided throughout the process. Program with 10 training modes available. Intelligent detection of whether the electrode pads are pasted. With a pause function, the volume can be adjusted.

The standard configuration:

An advanced resuscitation full-body manikin

One luxury hand-pulled human body hard plastic box

A computer (user-configured or optional)

One automatic external defibrillator (for training)

A set of CPR installation software
A resuscitation pad
Barrier mask (50 sheets/box) a box
Four sets of interchangeable lung sac devices
Two interchangeable skins
2010 operation guide CD-ROM



Samco Medical

Website: www.samcomedical.com,