



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Semi Fowler Hospital Bed

Product Code :
SAMCOMEDICAL1001A003

Description :

Semi Fowler Hospital Bed

Technical Specification :

Size : 198 L x 90 W cm x 60 H cms. Frame work made of strong rectangular steep tubes. Two section of CRCA M.S. Sheet uniformaly perforated. Backrest section maneuvered by a screw lever from foot end. Head and leg bows of steel tubular pipes assembled to laminated pannels. Legs of ERW tubular peipes mounted on four heavy duty rubber shoes. Pre-treated and epoxy powder coated. Optional Accessories : * Side Railing * Saline rod * Mosquito Net Poles Note : Above Bed can also supplied with castors.



Samco Medical

Website: www.samcomedical.com,