



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
EGG Incubator

Product Code :
SAMCOMEDICAL270001A001

Description :

EGG Incubator Backed by our efficient workforce, we are able to offer our clients with EGG Incubator. In our products, there is a gap between inner and outer wall that is filled with the superior quality mineral glass wool. EGG Incubator Temperature in these products can be reached from ambient to 70½ C, which is controlled by the ½ EGO½ or ½ JUMO½ German capillary type thermostat having an accuracy of ½ 0.5½ C. In addition to this, EGG Incubator our products have manual tilting device that is useful in tilting the tray holding the eggs. EGG Incubator Capacity: ½ 50 eggs ½ 100 eggs ½ 200 eggs ½ 250 eggs EGG Incubator, Lab Incubator Exporters, Lab Incubator Manufacturers



Samco Medical

Website: www.samcomedical.com,