



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hybridization Incubator

Product Code :
SAMCOMEDICAL270001A002



Description :

Hybridization Incubator Our company specializes in offering qualitative range of Hybridization Incubator. The inside portion of the product offered by our organization is made of stainless steel. Hybridization Incubator These products have the temperature which ranges above ambient $5\frac{1}{2}^{\circ}\text{C}$ to $90\frac{1}{2}^{\circ}\text{C}$. our product has a circulating fan ensuring the uniformity of temp. It also has a carrier to hold bottles which can be easily on front of rotator mechanism and can also be easily taken out. Hybridization Incubator bottles are 6 in number and of 60 mm Dia x 200 mm depth. Hybridization Incubator carrier rotates with bottles parallel to the door in a variable speed which helps in viewing the glass window in the specimen inside the chamber. Our product is complete with the microprocessor PID digital temp. indicator-cum-controller having Pt 100 sensor. Hybridization Incubator Features: Hybridization Incubator Leak proof bottles Hybridization Incubator Double walled Hybridization Incubator Inside stainless steel (SS) Hybridization Incubator With Viewing glass window Hybridization Incubator, Lab Incubator Exporters, Lab Incubator Manufacturers



Samco Medical

Website: www.samcomedical.com,