

**Product Name :**  
Walk In Incubator

**Product Code :**  
SAMCOMEDICAL270001A003



#### Description :

Walk In Incubator We hold expertise in undertaking a wide range of Walk In Incubator. These products are widely used for various incubation applications. Our products can be easily installed even in the hallway room or wherever convenient. Walk In Incubator are designed with the enameled mild steel sheet exterior, aluminum/stainless steel, interior ventilators, internal lighting, exhaust fan among others. Our products are suitable for the insulation of 300 mm sq. inspection window at the door which permits inside view without opening the door. Walk In Incubator products is complete with the panel and bare prefabricated at the factory assembled at site. It also has cord and plug to work on 220 V single phase 50 Hz power supply. Walk In Incubator temperature of these products range from 5°C above ambient to 60°C 1°C and are controlled by the digital indicator-cum-controller. These are of three types. Available size W x H x D: 1.2x2.3x1.2 meter 2.3x2.3x2.3 meter These are same as but with controlled humidity up to 95% 5% and are available in the following size W x H x D: 1.2x2. 3x1.0 meter 2.3x2.3x.3 meter These are same as but temperature range 10°C to 60°C and are available in the following sizes W x H x D: 1.2x2. 3x1.0 meter 2.3x2.3x.3 meter Walk In Incubator, Lab Incubator Exporters, Lab Incubator Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),