



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Rotary Senior Microtome

Product Code :
SAMCOMEDICAL280001A002



Description :

Rotary Senior Microtome We have established ourselves as a renowned manufacturer and supplier of Rotary Senior Microtome. All our products are widely acclaimed by the customers for their corrosion & abrasion resistance, sturdy construction, low maintenance and easy handling. These products are designed and developed using premium raw material, which is sourced from the reliable vendors of the market. Technical specifications: $\frac{1}{2}$ 1-50 microns $\frac{1}{2}$ Total feed excursion 28 mm with opening of object clamps 37x27 mm $\frac{1}{2}$ Overall dimensions 16 x 10 x 9 $\frac{1}{2}$ $\frac{1}{2}$ Complete with razor 120 mm with 120 mm with back & handle $\frac{1}{2}$ Honing plate $\frac{1}{2}$ Object holder set of three $\frac{1}{2}$ Oil can $\frac{1}{2}$ Packed in a wooden box , Microtome Instrument Exporters, Microtome Instrument Manufacturers



Samco Medical

Website: www.samcomedical.com,