



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Precision Rotary Microtome Lipshaw Type

Product Code :
SAMCOMEDICAL280001A004



Description :

Precision Rotary Microtome Lipshaw Type We are instrumental in offering a wide range of Precision Rotary Microtome Lipshaw Type. These products are designed with elegant and sturdy design using the latest technology. Our products are complete with counter balancing wheel, automatic feed mechanism and feed release. These products has light aluminum cover hinged at back which protects the the internal mechanism. Our products has a heavy casted case plate fitted with heavy rubber pads providing better grip. These products provides smooth operation for uniform serial sectioning and adjustable from 1 to 50 microns instep of 1 micron. feed mechanism is independent of other movements. These are supplied completely with one knife having back and handle in wooden box, 6 block holders, honing plate with abrasive powder, one bottle of lubricating oil and rexine cover. Features: $\frac{1}{2}$ Independent feed mechanism with automatic safety device $\frac{1}{2}$ Sturdy design $\frac{1}{2}$ Longer service life $\frac{1}{2}$ High performance Precision Rotary Microtome Lipshaw Type, Microtome Instrument Exporters, Microtome Instrument Manufacturers



Samco Medical

Website: www.samcomedical.com,