



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Microtome Cryostat Automatic

**Product Code :**  
SAMCOMEDICAL280001A005



**Description :**

Microtome Cryostat Automatic With the help of our experienced engineers, we are able to design and develop a wide range of Microtome Cryostat Automatic. These products provides a rapid and easy means for the sectioning of natural tissue from 2 to 20 microns in step of 2 microns without distortion for quality microscopic examination. These products have a higher resolution of cryptological details required for the pathological diagnosis, histochemistry, fluorescence, biochemistry, microscopy among others. These products are simple and easy to learn and withstand condensation inside the chamber. Our products are made of entirely rust-proof material fitted with anti roll device to avoid curling of sections and are highly adjustable. These works in the temperature ranging up to  $-30\pm\frac{1}{2}$  C as selected by imported thermostat and maintained by hermetically sealed compressor system. Features:  $\pm\frac{1}{2}$  Easy visibility  $\pm\frac{1}{2}$  Easy operation Control panel includes:  $\pm\frac{1}{2}$  Thermostat control  $\pm\frac{1}{2}$  Temperature indicator i.e. Dial thermometer automatic timer  $\pm\frac{1}{2}$  Pilot lamp  $\pm\frac{1}{2}$  Line switch  $\pm\frac{1}{2}$  Line indicator  $\pm\frac{1}{2}$  Voltmeter which indicates the incoming voltage Microtome Cryostat Automatic, Microtome Instrument Exporters, Microtome Instrument Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),