



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Automatic Microtome Knife Sharpener

Product Code :
SAMCOMEDICAL280001A008



Description :

Automatic Microtome Knife Sharpener Our company specializes in offering qualitative range of Automatic Microtome Knife Sharpener. These products have the unmatched and the most skilled hard honing and stropping. Our products permits to set the extract angle of bevel of knife edge. These products have a constant speed motors to rotate the ground glass plate 500 mm Dia x 12 mm thick at constant speed. Salient features of the unit include the following: $\frac{1}{2}$ Easy operation $\frac{1}{2}$ Automatic reversing $\frac{1}{2}$ Eliminates the need for resurfacing the glass plate as the knife execution $\frac{1}{2}$ Full diameter of the rotating plate $\frac{1}{2}$ Easy cleaning and maintenance $\frac{1}{2}$ Easily sharpened Optional: $\frac{1}{2}$ Micro Coarse & Fine Abrasive (set of two bottles of 200 gm each) $\frac{1}{2}$ Spare glass plate $\frac{1}{2}$ Spare knife holder Automatic Microtome Knife Sharpener, Microtome Instrument Exporters, Microtome Instrument Manufacturers



Samco Medical

Website: www.samcomedical.com,