



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Semi Automatic Rotary Microtome

Product Code :
SAMCOMEDICAL280001A010



Description :

Semi Automatic Rotary Microtome We are offering Semi Automatic Rotary Microtome catering to the needs and requirements of our esteemed customers. This equipment is featured with a standard drive system for smooth sectioning and trimming. We have also installed LCD screen displays section, various trimming thickness range and number of sections in this microtome. There is a lockable hand wheel, interchangeable specimen clamp along with a safety alarm system and huge volume waste tray for easy operation. Semi Automatic Rotary Microtome calibrated controls assure an accurate positioning to make re-cuts possible. Specifications: $\pm 1/2$ Section thickness setting range: 0.25-60 $\pm 1/2$ Trimming thickness setting Range: 0.25-60 $\pm 1/2$ Setting values: 0.25-1 $\pm 1/2$ m increment step 0.25 $\pm 1/2$ m $\pm 1/2$ 1-10 $\pm 1/2$ m increment step 1 $\pm 1/2$ m $\pm 1/2$ 10-20 $\pm 1/2$ m increment step 2 $\pm 1/2$ m $\pm 1/2$ 20-60 $\pm 1/2$ m increment step 5 $\pm 1/2$ m $\pm 1/2$ Minimum setting value: 0.25 $\pm 1/2$ m Semi Automatic Rotary Microtome, Microtome Instrument Exporters, Microtome Instrument Manufacturers



Samco Medical

Website: www.samcomedical.com,