



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Semi Automatic Cryostat Microtome

**Product Code :**  
SAMCOMEDICAL280001A012



**Description :**

Semi Automatic Cryostat Microtome Customers can avail Semi Automatic Cryostat Microtome from us. This machine is featured with pro-environment refrigerant, a rapid refrigeration mode, LED Screen to display data including sample travel, temperature and section thickness. Semi Automatic Cryostat Microtome is also featured with easy to operate control panel. There is also a closed drainage system with this machine to assure collection of condensed fluids in a container. Features:   
• Working area is easy to clean and disinfect   
• Easy to read LED screen   
• Ergonomic design   
• Specimen retraction function protects the specimen from blade injury   
• Trimming / sectioning button makes it easy to switch to the desired mode   
• Rapid refrigeration mode allows desired temperature drop in 60 mins   
Semi Automatic Cryostat Microtome, Microtome Instrument Exporters, Microtome Instrument Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),