



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Laminar Air Flow Cabinet

**Product Code :**  
SAMCOMEDICAL290001A001

#### Description :

Laminar Air Flow Cabinet Our clients can avail from us Laminar Air Flow Cabinet that is featured with double filtration of air through coarse pre-filter (up to 5 microns) and Hepa filter (down to 0.3 microns). Our cabinets assure the filtration of both particulate and biological contaminants. In addition to this, our expert quality controllers thoroughly check all the products on various parameters to ensure flawlessness. Features are given below: 1/2 Fluorescent lighting 1/2 AUE or equivalent motor 1/2 An acrylic block type manometer 1/2 Whole cabinet made of commercial board & laminated with sunmica 1/2 Whole cabinet made of stainless steel sheet duly smooth polished enhances the life of cabinet Specifications: Working Size 2'1/2x2'1/2x2'1/2 3'1/2x2'1/2x2'1/2 4'1/2x2'1/2x2'1/2 6'1/2x2'1/2x2'1/2 8'1/2x2'1/2x2'1/2 Size of HEPA Filters 2'1/2x2'1/2x6'1/2 3'1/2x2'1/2x6'1/2 4'1/2x2'1/2x6'1/2 3'1/2x2'1/2x6'1/2 4'x2'1/2x6'1/2 No. of HEPA Filters (Supplied with main unit) 1 1 1 2 2 No. Of Pre Filter 1 1 1 2 2 Illumination 1x20 W 1x20 W 2x20 W 2x40 W 4x40 W U.V. Germicidal Light 1x1'1/2 feet 1x1'1/2 feet 1x3 feet 2x3 feet 2x3 feet Laminar Air Flow Cabinet, Clean Air System Exporters, Clean Air System Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),