



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Biological Safety Cabinet

Product Code :
SAMCOMEDICAL290001A002



Description :

Biological Safety Cabinet In our large line of products we are offering Biological Safety Cabinet. This cabinet is fabricated using quality-approved thick board duly laminated sheet clad. We also use heavy gauge stainless steel sheet in this work table. Biological Safety Cabinet side panels of this tables is made of thick transparent flexi glass which is duly framed. It is featured with pre-filter and HEPA filter. Biological Safety Cabinet size of the HEPA filter is size 1'6" x 1'6" x 6'6". In addition to this, we check these cabinets on the basis of various well-defined parameters. Features: 1'6" U.V. germicidal Light 1'6" static pressure manometer 1'6" Exhaust system (suitable for 6 feet duct) 1'6" Virus burn out unit Biological Safety Cabinet, Clean Air System Exporters, Clean Air System Manufacturers



Samco Medical

Website: www.samcomedical.com,