



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Horizontal Laminar Air Flow Cabinet

**Product Code :**  
SAMCOMEDICAL290001A003



**Description :**

Horizontal Laminar Air Flow Cabinet We are providing Horizontal Laminar Air Flow Cabinet catering to needs of our esteemed customers. This cabinet is featured with table top model and the floor is not mounted. Our offered cabinet can be placed on existing laboratory platform. Featured with static pressure manometer, U.V. fluorescent light and cork for gas & air, this cabinet is highly appreciated by our customers. Optional: Whole cabinet made of commercial board & laminated with sunmica Whole cabinet made of stainless steel sheet duly smooth polished enhances the life of cabin Specification: Working Size 2x2x2 3x2x2 4x2x2 6x2x2 8x2x2 Size of HEPA Filters 2x2x6 3x2x6 4x2x6 3x2x6 4,x2x6 No. Of Pre Filter 1 1 1 2 2 Illumination 1x20 W 1x20 W 2x20 W 2x40 W 4x40 W U.V. Germicidal Light 1x1 feet 1x1 feet 1x3 feet 2x3 feet 2x3 feet Horizontal Laminar Air Flow Cabinet, Clean Air System Exporters, Clean Air System Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),