



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Horizontal Laminar Air Flow Cabinet

Product Code :
SAMCOMEDICAL290001A003



Description :

Horizontal Laminar Air Flow Cabinet We are providing Horizontal Laminar Air Flow Cabinet catering to needs of our esteemed customers. This cabinet is featured with table top model and the floor is not mounted. Our offered cabinet can be placed on exiting laboratory platform. Featured with static pressure manometer, U.V. fluorescent light and cork for gas & air, this cabinet is highly appreciated by our customers. Optional: 1½ Whole cabinet made of commercial board & laminated with sunmica 1½ Whole cabinet made of stainless steel sheet duly smooth polished enhances the life of cabin Specification: Working Size 21½x21½x21½ 31½x21½x21½ 41½x21½x21½ 61½x21½x21½ 81½x21½x21½ Size of HEPA Filters 21½x21½x61½ 31½x21½x61½ 41½x21½x61½ 31½x21½x61½ 4,x21½x61½ No. Of Pre Filter 1 1 2 2 Illumination 1x20 W 1x20 W 2x20 W 2x40 W 4x40 W U.V. Germicidal Light 1x11½ feet 1x11½ feet 1x3 feet 2x3 feet 2x3 feet Horizontal Laminar Air Flow Cabinet, Clean Air System Exporters, Clean Air System Manufacturers



Samco Medical

Website: www.samcomedical.com,