



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Inoculation Chamber

Product Code :
SAMCOMEDICAL290001A005



Description :

Inoculation Chamber Our numerous customers can avail Inoculation Chamber. This chamber is used for transferring tissues and other tissues culture work .It is used as dry-box. We use qualitative ply board clad in the manufacturing process of this cabinet. Inoculation Chamber inner side is exposed and this area is painted with chemically resistant epoxy. This cabinet is featured fluorescent lighting and germicidal ultra violet tube. Chamber size: 100x60x90 cm 100x60x60 cm 70x50x40 cm Chamber size (Inoculation chamber economy): 100x60x60 cm 70x50x40 cm Inoculation Chamber, Clean Air System Exporters, Clean Air System Manufacturers



Samco Medical

Website: www.samcomedical.com,