



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Grossing Table

**Product Code :**  
SAMCOMEDICAL290001A007



**Description :**

Grossing Table Customers can avail Grossing Table from us. The blowers of this table can suck the the contaminated air & formalin vapors. Grossing Table table is featured with lighting system, metric ruler, shelf, swiveled-magnifier, SS table, Sink & Tap, drawers, storage compartment under table, Towel / Tissue paper stand and waste bin. We also check these products on the basis of various parameters. Working Area (L $\times$ W $\times$ H): 4'  $\times$  3'  $\times$  6' 5'  $\times$  2.5'  $\times$  3' Optional Accessories: Built in U.V. germicidal light Exhaust duct upto 6 Grossing Table, Clean Air System Exporters, Clean Air System Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),