



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Fume Hood

Product Code :
SAMCOMEDICAL290001A008



Description :

Fume Hood The Fume Hood offered by us fabricated using wooden board. It is specially designed to throw-toxic and fumes and it protects costly instruments. Fume Hood outer side is finished with sunmica and inner side is epoxy painted. Fume Hood front facing panel of this product is fitted with 15/5 Amp socket with switches for exhaust system and fluorescent light. Working Chamber Size: 3' 1/2 x 2' 1/2 x 2' 1/2 4' 1/2 x 2' 1/2 x 2' 1/2 5' 1/2 x 2' 1/2 x 2' 1/2 6' 1/2 x 2' 1/2 x 2' 1/2 Optional: High CFM exhaust system (for duct exceeding 10 feet in length) S.S. lining on work bench. PVC fume duct (perft.) Fume Hood, Clean Air System Exporters, Clean Air System Manufacturers



Samco Medical

Website: www.samcomedical.com,