



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Resuscitation Unit

**Product Code :**  
SAMCOMEDICAL300001A001

**Description :**

Resuscitation Unit The controller is based on Micro computerized system i.e. Servo and Manual Modes i.e. Self test system provided i.e. Large bright displays are easy to monitor and for viewing from a distance i.e. Double digital display system provided to measure the Infant temperature and to set the required temperature i.e. Feather touch keys provided for easy operation i.e. Bar graph LED system provided to indicate the heater output i.e. Fault indication and alarms for i.e. Baby probe failure i.e. Low temperature i.e. High temperature cut off i.e. Heater failure i.e. Power failure i.e. Baby probe is made of highly sensitive sensor and interchangeable which does not need field calibration, the probe is made up of silicon sheathed cable, which allows easy cleaning i.e. Single element imported CE marked ceramic heater placed in a parabolic reflector i.e. The large baby tray made of stainless steel frame with three side acrylic collapsible walls i.e. IV stand provided i.e. Trendelenburg / fowler system provided i.e. Unit is powder coated i.e. Observation lamp provided i.e. Monitor tray provided Resuscitation Unit, Infant Incubator Exporters, Infant Incubator Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),