



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Suction Apparatus Oil Immersion Type

**Product Code :**  
SAMCOMEDICAL310001A002



**Description :**

Suction Apparatus Oil Immersion Type Innovation is our forte and for this, we provide our clients with Suction Apparatus Oil Immersion Type. All our products are widely acclaimed by the customers for their easy installation, low maintenance, excellent performance and longer service life. These products are designed and developed using high-grade raw material, which is procured from the reliable vendors of the market. Features: Suction Mechanism is oil immersed Suction Capacity is high Two 1500 CC bottles having Chrome Plated Metal Caps Automatic Cut off preventing overflow Stainless steel top Regular Valve to control speed of suction Made for running continuously for many hours Moves on quality casters wheels Elegant cabinet Foot Switch available as on optional extra Available also in stainless body Pilot lamp Noiseless Operation Filter is provided to absorb moisture and water particles entering into the rotor Additional Motor Protection Switch For operation on 230 Volts, 50 cycles, Single Phase, A.C Suction Apparatus Oil Immersion Type, Suction Unit Exporters, Suction Unit Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),