



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Medical Suction Units

**Product Code :**  
SAMCOMEDICAL310001A003



**Description :**

Medical Suction Units We are instrumental in offering our customers with Medical Suction Units that are designed in accordance with international quality norms & standards. These products are widely acclaimed for their excellent performance, low maintenance, easy installation, zero defect and longer service life. Apart from this, we provide these products in different models, sizes and designs at market leading prices. Standard consumables include the following: Two autoclavable polycarbonate jars of 2000ml with overflow valve system One antibacterial and hydrophobic filter One complete set of silicone / PVC tubes One conical connector for Cannula One Power Cord Options: Available with 1.5 ltr., 2 ltrs. Glass Jars Technical specification: Motor: Oil less and maintenance free piston pump Power Feeding: 220V  $\pm$  10% 50Hz Maximum Vacuum- (Adjustable): 700mmHg  $\pm$  10mmHg Max free air flow rate: 30l/min Duty Cycle: Non-stop operating Years of warranty: 1 Shipping Carton: 1 Medical Suction Units, Suction Unit Exporters, Suction Unit Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),