



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Suction Machine

Product Code :
SAMCOMEDICAL310001A004



Description :

Suction Machine ½ One autoclavable polycarbonate jar with overflow valve system ½ One antibacterial and hydrophobic filter ½ One complete set of silicone tubes ½ One conical connector for cannula ct ½ One Power Cord Technical specifications: ½ Motor: Oil less and maintenance free piston pump ½ Power Feeding: 220 50Hz other configurations on request ½ Maximum Vacuum (adjustable): -0.80bar -80kPa -600mmHg ½ Max free air Flow rate: 30l/min ½ Years of warranty: 1 ½ Shipping Carton: 4 Suction Machine, Suction Unit Exporters, Suction Unit Manufacturers



Samco Medical

Website: www.samcomedical.com,