



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Standard Consumables Suction Unit

Product Code :
SAMCOMEDICAL310001A005



Description :

Standard Consumables Suction Unit Two autoclavable polycarbonate jars of 2000ml with overflow valve system
1/2 One antibacterial and hydrophobic filter 1/2 One complete set of silicone / PVC tubes 1/2 One conical connector for cannula 1/2 One Power Cord Options: Available with Oil Free Piston Pump Technical specifications: 1/2 Motor: Rotary Pump 1/2 Power Feeding: 220V 50Hz 1/2 Other configurations on request
Maximum Vacuum (Adjustable): 700mmHg 1/2 Max free air flow rate: 50l/min 1/2 Duty Cycle: Non - stop operating 1/2 Years of warranty: 2 1/2 Shipping Carton: 1 Standard Consumables Suction Unit, Suction Unit Exporters, Suction Unit Manufacturers



Samco Medical

Website: www.samcomedical.com,