



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Portable Suction Machine

Product Code :
SAMCOMEDICAL310001A006



Description :

Portable Suction Machine Two autoclavable polycarbonate jars with overflow valve system ½ One antibacterial and hydrophobic filter ½ One complete set of silicone / PVC tubes ½ One conical connector for cannula ½ One Power Cord Options: Also available with Rotary Oil Dip Pump Technical specifications: ½ Motor: Oil less and maintenance free piston pump ½ Power Feeding: 220 50Hz other configurations on request ½ Maximum Vacuum (adjustable): -700mmHg+/-10mmHg ½ Max free air flow rate: 30l/min ½ Years of warranty: 1 ½ Shipping Carton: 1 Portable Suction Machine, Suction Unit Exporters, Suction Unit Manufacturers



Samco Medical

Website: www.samcomedical.com,