



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Digital Ultrasonic Cleaner

Product Code :
SAMCOMEDICAL320001A002



Description :

Digital Ultrasonic Cleaner Digital Ultrasonic Cleaner Innovation is our forte and we strive to design and develop a precision engineered range of Digital Ultrasonic Cleaner. We design and develop all the products using high-grade raw material and latest technology to ensure that final output meets with international quality standards & norms. In addition to this, we provide customization as per the specifications laid down by the customers.

Specifications: Tank Dim.: 45cm x 15cm x 15cm Liter Cap.: 10 Ltrs U/S Power: 200 Watts U/S Freq.: 25 khz Input:180 Vac - 270 Vac Warranty: 1 Year S.M.P.S. ensures smooth working of cleaners within the voltage variation of 180-270 vac. Half the weight & size of the conventional machines Competitive rates efficient service backup of India Models can be made as per your specific requirements Heater ensures high efficiency in cleaning Cleaners are equipped with Imported Transducers for smooth and long life working. Industries Applicable:- Hospitals, Nursing Homes, Educational Institutes, Electronic Industry, Jewelry Industry, Hotel Industry, Pharmaceutical Industry & Automobile Industry Digital Ultrasonic Cleaner, Ultrasonic Cleaner Exporters, Ultrasonic Cleaner Manufacturers



Samco Medical

Website: www.samcomedical.com,