



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Cold Light Sources

Product Code :
SAMCOMEDICAL330001A001

Description :

Cold Light Sources We are engaged in providing our clients with Cold Light Sources that are available in different models and sizes at market leading prices. These products are widely appreciated by the customers for their corrosion & abrasion resistance, excellent performance, low maintenance and longer functional life. Furthermore, our products can be customized as per the specifications laid down by the clients. Technical specifications: 1/2 Model: Single lamp/double lamp 1/2 Output power: 250 watts (24 volts) 1/2 Input: 180 Vac - 270 Vac Features: 1/2 S.M.P.S. based design ensures smooth working of light source within the voltage variation of 180 -270 vac 1/2 Half the weight & size of the conventional light sources. 1/2 Light intensity controlled with smooth rotary potentiometer Cold Light Source, Cold Light Source Exporters, Cold Light Source Manufacturers



Samco Medical

Website: www.samcomedical.com,