

Product Name :
Medical Refrigerator

Product Code :
SAMCOMEDICAL350001A002



Description :

Medical Refrigerator Medical/Blood Bank Refrigerators high efficiency compressor for silent working and energy saving. ½ Medical/Blood Bank Refrigerators Highly effective condenser and expansile evaporator provide quick freezing. ½ "ATICO" Low Cost Medical Refrigerator is a specialized refrigerating equipment for cold storage of pharmaceuticals as well as biological products and so forth. ½ It is suitable for hospitals, drugstore, pharmaceuticals factories, sanitation and anti epidemic stations and clinics. ½ Microprocessor based temperature control system; combined with two exact sensors to defrost. Automatically turn on the compensatory heater based on the temperature change. A stable and accurate chamber temperature can be maintained at ambient temperature ranging from 0½C to 32½C. ½ Forced air circulation system with intelligent control ensures excellent temperature uniformity. ½ Automatic temperature control : digital display of temperature. The temperature from 2 to 10½C can be set freely, controlling precision 1½C, display accuracy 0.1½C ½ Perfect alarm systems, which provider numerous audible/visual alarm functions in case of high or low temperature condition, sensor abnormality, door open & power failure to ensure the safety of storage. ½ Exterior made from high quality structural steel plate, painted by advanced coating technology after the anti corrosion treatment of phosphorization, the superficial luster is gentle, interior made from ABS engineering plastics. ½ Double deck glass door, inert gas, equipped with special lock, four casters are mounted under the bottom. ½ Interior fluorescent lighting provides a clear view of stored items. ½ Shelves made of quality steel wire soaked with plastic provide convenient access to the stored articles and easy cleaning. ½ High efficiency compressor for silent working and energy saving. ½ Highly effective condenser and expansile evaporator provide quick freezing Medical Refrigerator, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers

Samco Medical

Website: www.samcomedical.com,

