



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Blood Collection Monitor

**Product Code :**  
SAMCOMEDICAL350001A005



**Description :**

Blood Collection Monitor Volume can be set 1 ml increment  
• Facility of pausing during collection and also can change preset volume  
• Micro controller based program to give better accuracy  
• Auto tare facility  
• Motorized pinching system at the end of collection to prevent excess blood from the desired volume  
• Facility to check, set volume, gram and time totalizer during collection  
• Audio video alarm if the flow is less than 15 ml/ min and Audio Alarm at the end of collection  
Blood Collection Monitor, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),