



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Blood Bag Tube Sealer

Product Code :
SAMCOMEDICAL350001A006



Description :

Blood Bag Tube Sealer Oscillation Frequency 40.68 Mhz Radio Frequency output 15W Sealing time Approx. 1.5 sec. Input Volt 230V AC Power Frequency 50 Hz Power Consumption 100 VA Weight 5.5 Kg approx Size 180 mm x 215 mm x 245 mm (HxDxW) Blood Bag Tube Sealer, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers



Samco Medical

Website: www.samcomedical.com,