



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Platelet Incubator

**Product Code :**  
SAMCOMEDICAL350001A007



**Description :**

Platelet Incubator Micro Controlled Temperature controller with Audio Visual Alarm for high & low temperature  
½ Display of set and process value ½ Extra - strength non CFC refrigeration system ½ The unit is mounted on castor wheel for easy mobility ½ CFC Free PUF (Poly Urethane Foam) Insulation to ensure temperature stability and save power consumption Optional Accessories: ½ 7 Days circular temperature Recorder ½ Multi Channel Data logging system Platelet Incubator, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),