



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Plasma Extractor

**Product Code :**  
SAMCOMEDICAL350001A011



**Description :**

Plasma Extractor Plasma Extractor is designed to Extract Blood component from centrifuged bags ½ The spring loaded front panel applies pressure on the collected bag causing the liquid to contain in transfer bag ½ Its sturdy construction makes it durable and reliable, easy to sue and portable ½ The unit is available in both versions automatic and manual Technical Specification: Sensor - Infra Red Clamping Control Manual Motor Activated Pressure Plate Acrylic Alarm - Audio Visual Plasma Extractor, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),