



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Blood Counting Chamber

Product Code :
SAMCOMEDICAL360001A002

Description :

Blood Counting Chamber have counting grids which are engraved into the chamber bottom, under the microscope the counting grid Counting chambers are instruments for precise measuring.

Technical Specification :

Blood Counting chambers is made of special optical glass and serve to determine the number of particles in liquids, standard counting chambers have counting grids which are engraved into the chamber bottom, under the microscope the counting grid appears as a grid of dark lines. Blood Counting chambers have counting grids which are engraved into an annealed metal coating, under normal microscope settings the counting grid appears as a grid of bright lines, by altering the contrast, the counting grid can also be made to appear dark. Blood Counting Chamber China, Medical Educational Aids Exporters China, Medical Educational Aids China Manufacturers



Samco Medical

Website: www.samcomedical.com,