



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Aneroid Adult Size Blood Pressure

Product Code :
SAMCOMEDICAL360001A003

Description :

Aneroid Sphygmomanometers, For Hospital, Clinic Training. Aneroid Adult Size Blood Pressure

Technical Specification :

Aneroid Sphygmomanometers, For Hospital, Clinic Training

Professional quality Aneroid Sphygmomanometer Blood Pressure Cuff. Features a standard inflation bulb, standard air release valve and calibrated navy nylon cuff with velcro closure, to insure a secure fit and accurate reading. Applies the mechanical manometer to measure the blood pressure Easy to maintain Works equally well in hospitals, nursing homes or home health care. Medical Educational Training Aids Exporters, Medical Educational Training Aids Manufacturers



Samco Medical

Website: www.samcomedical.com,