



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ESR Stand Stainless Steel

Product Code :
SAMCOMEDICAL360001A004

Description :

ESR Stand Stainless Steel, ESR Stand is push type which can fit 12 Tubes

Technical Specification :

ESR means Erythrocyte Sedimentation Rate.

ESR Stand is also known as Westergren Stand, ESR Westergren Stand and ESR Westergren Apparatus.

ESR Stand is push type which can fit 12 Tubes.

ESR Stand 12 Tubes is made of Stainless Steel.

ESR Stand is a simple non-specific screening test that indirectly measures the presence of inflammation in the body associated with conditions such as infections, cancers, and autoimmune diseases.

ESR reflects the tendency of red blood cells to settle more rapidly in the face of some disease states because of increases in plasma fibrinogen, immunoglobulins, and other acute-phase reaction proteins.



Samco Medical

Website: www.samcomedical.com,