



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Vertical Foot Suction

Product Code :
SAMCOMEDICAL390001A001

Description :

Vertical Foot Suction We are highly appreciated in the international market for our exquisite collection of Vertical Foot Suction. These products are manufactured using superior raw material, which is procured from the reliable vendors of the market. In addition to this, our products are stringently checked by expert quality controllers to ensure flawlessness. Other information: ½ Portable and non string pedal operated unit ½ One wide mouth glass jar of approx. 1600 ml ½ Piston type pump is very effective to create the vacuum instantly ½ Out right vacuum 600 mm Hg ½ Self sealing rubber lids Vertical Foot Suction, Suction Apparatus Exporters, Suction Apparatus Manufacturers



Samco Medical

Website: www.samcomedical.com,