



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Foot Operated Electric Suction Apparatus

**Product Code :**  
SAMCOMEDICAL390001A002



**Description :**

Foot Operated Electric Suction Apparatus Our organization is counted amongst the most distinguished manufacturers and exporter of Foot Operated Electric Suction Apparatus. All our products are widely acclaimed by the customers for their corrosion & abrasion resistance, dimensional accuracy, sturdy construction and low maintenance. Also, our products are available in different sizes, designs and models at market leading prices. Other information: ½ Our foot suction apparatus is most suitable for places where electric supply is uncertain ½ Well-designed cabinet made of mild steel epoxy powder coated with stainless steel top provided ½ Two separate Suction System electric suction apparatus and vertical foot operated suction ½ 1 / 4 HP heavy duty motor fitted to oil immersed vacuum pump ½ Two wide mouth glass bottles of approx. 1600 ML each ½ Synthetic Airtight Rubber Lids. Non-collapsible PVC Tubing ½ Two separate vacuum gauge for each system ½ Vacuum of electric suction 700 mm + 10mm Hg in our electric suction apparatus ½ Ultimate vacuum of vertical foot suction 600 mm Hg Foot Operated Electric Suction Apparatus, Suction Apparatus Exporters, Suction Apparatus Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),