



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Suction Apparatus Trolley Model

Product Code :
SAMCOMEDICAL390001A003



Description :

Suction Apparatus Trolley Model Ideal & suitable for Hospitals, Nursing Homes and Clinics
• Cabinet made of mild steel sheet with epoxy powder coated
• 1/2 HP heavy duty motor. Oil immersed vacuum pump
• Two wide mouth glass bottles of approx. 1600 ml. each
• Synthetic air tight rubber lids. Overflow cut-off device
• Non-collapsible PVC Tubing
• Vacuum-700 +/- 10mm Hg
• Trolley mounted on friction free wheels
Suction Apparatus Trolley Model, Suction Apparatus Exporters, Suction Apparatus Manufacturers



Samco Medical

Website: www.samcomedical.com,