



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dual Bottle Suction Machine

Product Code :
SAMCOMEDICAL390001A004



Description :

Dual Bottle Suction Machine Each suction method independent of other
• Comes with separate vacuum gauge support with reading from 0 to -760 mmHg
• Electric unit option coming fitted with suitable motor support that provides for continuous flow rate
• Unit comes with provision of foot operation in case of power failure
• Bottle coming fitted with automatic float cut-out so as to protect pump from accidentally sucking in fluids
• Front panel comes with support of pilot indicator lamp, ON/OFF switch, vacuum gauge and vacuum control knob
Dual Bottle Suction Machine, Suction Apparatus Exporters, Suction Apparatus Manufacturers



Samco Medical

Website: www.samcomedical.com,