



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Horizontal Cylindrical Sterilizer

Product Code :
SAMCOMEDICAL400001A002



Description :

Horizontal Cylindrical Sterilizer Steam Working Pressure 1.26 kgf/cm² (2.2 kgf/cm² in case High speed sterilizer
Steam Working Temperature 121oC Operating Voltage 240, single phase, AC supply, 50 Hz ½ Sterilization
Period: 45 to 50 minutes ½ Radiant locking system, cold water drainage system ½ Steam Exhaust: 5 to 7
provided with timer from 1 to 60 Min with double safety valve, steam releasing valve, pressure gauge for the
pressure inside the chamber, and one pressure gauge to show the pressure of outer chamber both the chamber
made of stainless steel lid of the unit made of heavy duty gun metal piece for extra safety, gasket, water level
indicator to see the water inside, heavy duty power plug with socket provided with the unit, accessories: one
spare element, and two gasket , our common size 400x1200 mm Horizontal Cylindrical Sterilizer, Medical
Sterilizers Exporters, Medical Sterilizers Manufacturers



Samco Medical

Website: www.samcomedical.com,