



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Radiant Warmer

Product Code :
SAMCOMEDICAL410001A001

Description :

Infant Radiant Warmer We are one of the preeminent organizations, involved in manufacturing and exporting a wide range of Infant Radiant Warmer. These products are widely acclaimed by the customers for their longer service life, excellent performance, low maintenance and easy handling. Furthermore, we have the capability to customize all the products as per the specifications laid down by the customers. Control Panel : Infant Warmers are designed to provide advance warming capabilities while using one touch operations to simplify procedures. Advance microprocessor controlled system with self check on the electronic circuitry skin / air / manual control. Digital display of baby temperature / air temperature . Soft touch display of control temperature. Manual / servo two function control modes provide the clinician more flexibility while optimizing the thermal environment. In addition the skin / Air temperature is larger than the other set temperature display. Skin temperature can be converted from degree C to degree F by the push of a button. Baby Bed : Fixed The spacious baby bed incorporates easily drop down sides for maximum access to the patient. Soft cushion mattress for the comfort of the little patient. Alarms : The control panel is also the heart of the system for a completely integrated alarm system with both audible and visual indicators in the event of Probe fail Skin high Skin low Power fail Heater fail Manual Mode : Facility to converting the warmer from servo mode to manual mode for non-stop working. In Manual mode warmer will start from 40% heat and can be increased upto 99% according to usage. Heater Assembly : The heater assembly located above the center of the support structure consists of radiant heater and reflector with inbuilt halogen examination lamp for assessment of infant color Infant Radiant Warmer, Infant Care Products Exporters, Infant Care Products Manufacturers



Samco Medical

Website: www.samcomedical.com,