



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Baby Incubator

Product Code :
SAMCOMEDICAL410001A005



Description :

Infant Baby Incubator The temp controlled unit is operated by advanced micro computer technology. The three digital display indication provide the Set Temperature, Air Temperature and Skin Temperature. The safety features are self test and alarm system. Which will detect the different failures. The removable heater and fan unit make it easy to clean and service. In addition, the Skin/Air temperature display is larger than the other set temperature display. Full visibility & accessibility of the baby is provided by the designed of acrylic hood. A full complement of alarm provides safety convenience, with both audible and visual indicator.

Electrical Specification: Power 200V/50Hz Power Consumption 300Watts Chassis Leakage >500 Amp Current Fuse 5Amp

Infant Baby Incubator, Infant Care Products Exporters, Infant Care Products Manufacturers



Samco Medical

Website: www.samcomedical.com,