



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
CFL Phototherapy Unit

Product Code :
SAMCOMEDICAL410001A006



Description :

CFL Phototherapy Unit CFL Provides the light intensity of 18 pw/cm2/nm (approx) at 45cm from Baby cot
Choice of 6 Blue or 4 Blue + 2 white CFL lamps for treatment and inspection
Height 120 cm to 160 cm from the base
Head Tilting Facility from -1800 + 1800 (Continuous)
3 swivel casters fixed on sleek T-shaped base
Powder coated scratch and rust proof body
Special mirror Coated Reflector For better reflect CFL Phototherapy Unit, Infant Care Products Exporters, Infant Care Products Manufacturers



Samco Medical

Website: www.samcomedical.com,