



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Double Surface Phototherapy Unit

**Product Code :**  
SAMCOMEDICAL410001A007



**Description :**

Double Surface Phototherapy Unit i.e. ½ Baby Bassinet with fold down panels i.e. ½ Head down position, on both sides i.e. ½ Phototherapy Unit fitted at a distance of 300 mm from baby bed i.e. ½ 4 Blue + 2 white tubes in Overhead Phototherapy Unit Choice of two kinds of tubes for use with i.e. ½ 4 nos. of blue light tubes with higher intensity of light in therapeutic range and 2 nos of white light tubes i.e. ½ 6 nos. of white light tubes with the facility to use 4 or 6 tubes at a time i.e. ½ Compact & mobile



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),