



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Single Surface Phototherapy Unit

Product Code :
SAMCOMEDICAL410001A012



Description :

Single Surface Phototherapy Unit
• Overhead Photo therapy Unit
• Choice of two kinds of tubes for use with
• 4 nos. of blue light tubes with higher intensity of light in therapeutic range and 2 nos. of white light tubes
• 6 nos. of white light tubes with the facility to use 4 or 6 tubes at a time
• Facility of height adjustment, so as to fit over an incubator or baby cot etc
• Fold up facility during inspection of neonates or storage
• Angle folding facility, for using along with radiant warmer
• Compact & mobile Single Surface Phototherapy Unit,
Infant Care Products Exporters, Infant Care Products Manufacturers



Samco Medical

Website: www.samcomedical.com,