



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Harpenden Skinfold Calipers

Product Code :
SAMCOMEDICAL420001A001

Description :

Harpenden Skinfold Caliper The Harpenden Skinfold Caliper is a precision instrument designed for use in the performance of Skinfold thickness measurements (from which estimates of body fat are derived). The use of this instrument has been well established and documented over the past 40 years. Designed in 1958 in collaboration with D. J. M. Tanner, a prominent force in the use of Skinfold measurements in the derivation of body fat measurement, The Harpenden Skinfold Caliper has been used all over the world on projects of international importance and is the most prestigious Caliper available. It is the only Caliper CE marked under the Medical Devices Directive 93/42/EEC for a Class 1 Device with Measuring Function and is calibrated using masters traceable to National Standards. Dial Graduation: 0,20 mm Measuring Range: 0 mm to 80 mm Measuring Pressure: 10 gms/mm² (constant over range) Accuracy: 99.00% Repeatability: 0,20 mm Harpenden Skinfold Caliper, Harpenden Skinfold Caliper Exporters, Harpenden Skinfold Caliper Manufacturers



Samco Medical

Website: www.samcomedical.com,